


# Heuron AD

Heuron AD vs. FreeSurfer

## A. Accurate Segmentation

*Heuron AD uses a deep learning-based algorithm to accurately segment brain regions using T1-weighted MR images*

*Randomly selected, unbiased sample of MR database used to test the difference with **FreeSurfer***



Deidentified, anonymized MR data (Case A) is available for further assessment

# Heuron AD

Heuron AD vs. FreeSurfer

## B. Less Artifact


Heuron AD



FreeSurfer



Clearer visualization of  
cortical thickness and brain  
segmentation with  
**less artifact**




*Artifacts can lead to a misinterpretation and an error in volumetric assessment.*

# Heuron AD

Heuron AD vs. FreeSurfer

## C. 99 Segmented Brain Regions



Heuron AD

FreeSurfer

Midbrain, pons, and medulla can be analyzed separately

Only the brainstem can be analyzed

#No.	Label Name:	R	G	B	A
0	Unknown	0	0	0	0
2	Left-Cerebral-White-Matter	245	245	245	0
4	Left-Lateral-Ventricle	120	18	134	0
5	Left-Inf-Lat-Vent	196	58	250	0
7	Left-Cerebellum-White-Matter	220	248	164	0
8	Left-Cerebellum-Cortex	230	148	34	0
10	Left-Thalamus-Proper	0	118	14	0
11	Left-Caudate	122	186	220	0
12	Left-Putamen	236	13	176	0
13	Left-Pallidum	12	48	255	0
14	3rd-Ventricle	204	182	142	0
15	4th-Ventricle	42	204	164	0
17	Left-Hippocampus	220	216	20	0
18	Left-Amygdala	103	255	255	0
24	CSF	60	60	60	0
26	Left-Accumbens-area	255	165	0	0
28	Left-VentralDC	165	42	42	0
31	Left-choroid-plexus	0	200	200	0
41	Right-Cerebral-White-Matter	245	245	245	0
43	Right-Lateral-Ventricle	120	18	134	0
44	Right-Inf-Lat-Vent	196	58	250	0
46	Right-Cerebellum-White-Matter	220	248	164	0
47	Right-Cerebellum-Cortex	230	148	34	0
49	Right-Thalamus-Proper	0	118	14	0
50	Right-Caudate	122	186	220	0
51	Right-Putamen	236	13	176	0
52	Right-Pallidum	13	48	255	0
53	Right-Hippocampus	220	216	20	0
54	Right-Amygdala	103	255	255	0
58	Right-Accumbens-area	255	165	0	0
60	Right-VentralDC	165	42	42	0
63	Right-choroid-plexus	0	200	221	0
77	WM-hypointensities	200	70	255	0
173	Midbrain	242	104	76	0

⋮

99 Segmented Brain Structures

Region	HCP_vol	Heuron AD_vol
Left-Lateral-Ventricle	27925.6	26283.5
Left-Inf-Lat-Vent	1288.6	1053.9
Left-Cerebellum-White-Matter	10906.7	11427.7
Left-Cerebellum-Cortex	43225.4	43679.6
Left-Thalamus-Proper	5070.9	5267.6
Left-Caudate	2671.1	3028.7
Left-Putamen	2966.2	2881.6
Left-Pallidum	1283.3	1236.2
3rd-Ventricle	2206	1843.5
4th-Ventricle	1804.7	1559.1
Brain-Stem	18996.4	
Midbrain		5079.7
Pons		13350
Medulla		4157.9
Left-Hippocampus	2029.9	2235.4
Left-Amygdala	699.6	868.8
CSF	1230.8	1054.5
Left-Accumbens-area	267.3	209.5
Left-VentralDC	3036.6	1665.7
Left-vessel	13.2	
Left-choroid-plexus	883.9	1104.8
Right-Lateral-Ventricle	19651.5	18531.1
Right-Inf-Lat-Vent	847.5	784.4
Right-Cerebellum-White-Matter	10523.9	10973.4
Right-Cerebellum-Cortex	43713.9	44476.7
Right-Thalamus-Proper	5184.5	5129.7
Right-Caudate	2571	2780.8
Right-Putamen	3336.8	3266.1

⋮


Regional Volume Comparison  
FreeSurfer vs. Heuron AD

# Heuron AD

Heuron AD vs. FreeSurfer

## D. Cortical Thickness Measurement

Raw  
MR Data



Heuron AD

Heuron AD automatically assesses and calculates the cortical thickness information.

Recently, **cortical thickness** has been recognized  
as **valuable information for clinicians**



Alzheimer's & Dementia: Diagnosis, Assessment & Disease Monitoring ■ (2018) 1-10

Alzheimer's  
&  
Dementia

SPECIAL SECTION: State of the Field: Advances in Neuroimaging from the 2017 Alzheimer's Imaging Consortium

The personalized Alzheimer's disease cortical thickness index predicts likely pathology and clinical progression in mild cognitive impairment

Annie M. Racine<sup>a,b,c</sup>, Michael Brickhouse<sup>c</sup>, David A. Wolk<sup>c</sup>, Bradford C. Dickerson<sup>b,c,d,e,g</sup>,  
For the Alzheimer's Disease Neuroimaging Initiative<sup>1</sup>

<sup>a</sup>Aging Brain Center, Institute for Aging Research, Hebrew SeniorLife, Boston, MA, USA  
<sup>b</sup>Harvard Medical School, Boston, MA, USA  
<sup>c</sup>Frontotemporal Disorders Unit  
<sup>d</sup>Massachusetts Alzheimer's Disease Research Center, Department of Neurology, and Athinoula A. Martinos Center for Biomedical Imaging, Massachusetts General Hospital, Boston, MA, USA  
<sup>e</sup>Department of Neurology, Perelman School of Medicine, and Penn Memory Center, University of Pennsylvania, Philadelphia, PA, USA

**Abstract**  
**Introduction:** An Alzheimer's disease (AD) biomarker adjusted for age-related brain changes should improve specificity for AD-related pathological burden.  
**Methods:** We calculated a brain-age-adjusted "personalized AD cortical thickness index" (pADI) in mild cognitive impairment patients from Alzheimer's Disease Neuroimaging Initiative. We performed receiver operating characteristic analysis for discrimination between patients with and without cerebrospinal fluid evidence of AD and logistic regression in an independent sample to determine if a dichotomized pADI predicted conversion to AD dementia.  
**Results:** Receiver operating characteristic area under the curve was 0.69 and 0.72 in the two samples. Three empirical methods identified the same cut-point for pADI in the discovery sample. In the validation sample, 83% of pADI + mild cognitive impairment patients were cerebrospinal fluid AD biomarker positive. pADI + mild cognitive impairment patients (n = 63, 38%) were more likely to progress to AD dementia after 1 (odds ratio = 2.9) and 3 (odds ratio = 2.6) years.  
**Discussion:** The pADI is a personalized, magnetic resonance imaging-derived AD biomarker that predicts progression to dementia.  
© 2018 The Authors. Published by Elsevier Inc. on behalf of the Alzheimer's Association. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

**Keywords:** AD index; AD signature; Alzheimer's disease; Cortical thickness; Mild cognitive impairment

**1. Background**  
Positron emission tomography (PET) and cerebrospinal fluid (CSF) biomarkers are the gold standard for identifying individuals with molecular evidence of Alzheimer's disease (AD) neuropathology, but these procedures are invasive (CSF), expensive (PET), and only accessible in specialized centers (PET) [1,2]. Magnetic resonance imaging (MRI), on the other hand, is noninvasive, less expensive, and more readily available than PET but less specific than amyloid PET or CSF to AD-related neurodegeneration. Although the magnitude of hippocampal atrophy in patients scanned in vivo and followed to autopsy correlates with the burden of neurofibrillary tangle pathology [3], hippocampal atrophy can also be seen in patients with a variety of neurodegenerative and other pathologies [4–6]. Spatial patterns of regional

<https://doi.org/10.1016/j.dadm.2018.02.007>  
2352-8729/© 2018 The Authors. Published by Elsevier Inc. on behalf of the Alzheimer's Association. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

FLA 5.5.0 DTD ■ DADM239 proof ■ 17 March 2018 ■ 5:23 pm ■ ce

<https://www.sciencedirect.com/science/article/pii/S2352872918300150>